



Amalgam Separator Maintenance Log to be completed by Dental Staff Year: _____

Amalgam separators should be inspected weekly for drips or leaks. Any maintenance needed or performed should be noted. If the separator sludge collector or cartridge is made of clear plastic, note the level of sludge accumulation. Note the date you replace the collector box or cartridge and follow manufacturer instructions for recycling and obtaining a replacement collector box or cartridge.

Inspection Date	Initials	Leaks/Problems?	Maintenance Performed/ Comments

Name of Distributor of Amalgam Separator: _____
Contact Person: _____
Phone Number: _____